



RESUSCITATIVE TEE
PROJECT

8 HOUR DEEP DIVE

PRINCIPLES OF TEE OPERATION

CARDIAC ANATOMY

1:1 GOAL DIRECTED HANDS-ON PRACTICE

PROGRAM IMPLEMENTATION & LOGISTICS

PRACTICAL COMPETENCY ASSESSMENT

***“A STATE-OF-THE-ART COURSE
ON RESUSCITATIVE TEE”***

RESUSCITATIVETEE.COM



THE RESUSCITATIVE TEE WORKSHOP

PRESENTED BY
THE RESUSCITATIVE TEE PROJECT

Powered by



HEARTWORKS

RESUSCITATIVE TEE VIEWS		CARDIAC STRUCTURES	CLINICAL APPLICATION
ME4C 0-10 °		FOUR CHAMBERS, MV, TV PERICARDIUM	PERICARDIUM LV/RV FUNCTION WALL MOTION ABN VALVE LESIONS & DYSFUNCTION
MELAX 120-140°		LV, LA, MV AV	QUALITY CPR (AMC) LV FUNCTION PERICARDIUM MV/AV
TG MID SAX 0-20°		LV, RV	LV FUNCTION WALL MOTION ABN PERICARDIUM
ME BICAVAL 90-110°		IVC, RA, SVC	PROCEDURE GUIDANCE VENOUS GUIDEWIRE ECMO VOLUME RESPONSIVENESS
DESCENDING AORTA SAX 0°		DESCENDING AORTA	ARTERIAL GUIDEWIRE ECMO
DESCENDING AORTA LAX 90°		DESCENDING AORTA	ARTERIAL GUIDEWIRE ECMO

Design by @FTeranMD

GRAPHICS ©

